

SensiLase PAD-IQ®

The Intelligent Vascular Test!

	PAD-IQ	PAD3000
--	--------	---------

The SensiLase PAD-IQ is the ideal upgrade to the now discontinued PAD3000. The PAD-IQ offers a number of significant improvements over the PAD3000.



Skin Perfusion Pressure (SPP)	✓	✓
Pulse Volume Recording (PVR)	✓	✓
Ankle Brachial Index (ABI), assisted	✓	✗
SPP Test Sites / Exam Session	Multiple per Site; Up to 4 per Side; 8 per Session	Multiple per Site; Up to 4 per Side; 8 per Session
SPP Tests allows user to select low pressure inflation	✓	✗
PVR Test Sites / Exam Session	1 per Site; Up to 4 per Side; 8 per Session	1 per Site; Up to 4 per Side; 8 per Session
PVR test less sensitive to cuff wrap tightness	✓	✗
Simultaneous Bilateral Testing	✓	✗
Time to Perform Single Level Bilateral Test	7-10 minutes	20 minutes
Time to Perform Additional Bilateral Tests (per site)	4-6 minutes	10-15 minutes

©2017 VASAMED, Inc. All rights reserved. SensiLase PAD-IQ® and VASAMED® are registered trademarks of VASAMED, Inc.

VASAMED, Inc. 7615 Golden Triangle Drive - Suite A, Eden Prairie, MN 55345 Phone: 800.695.2737 Fax: 952.944.6022
04-00207-46-34 Rev 01




SensiLase PAD-IQ®

The Intelligent Vascular Test!

PAD-IQ PAD3000

Graphical User Interface



Angiosome-Based Test Site Selection	✓	✗
Intuitive Touch Screen User Interface	✓	✗
Supports Upper Extremity Documentation	✓	✗
Supports Wound Imaging / Documentation	✓	✗
Supports ICD-10 Codes	✓	✗
Portable All-In-One Design with Carrying Handle	✓	✗
Battery Life (AC Cordless)	4-6 hrs.	✗
Operating System	Win 7	Win XP
DICOM Store	✓	✗
Bayonet Air Hose Connecters	✓	✗

Contact us at 800-695-2737 or customerservice@vasamed.com to get upgraded to a new PAD-IQ today!